



Follow the table to determine the correct size.
Measurements are made 4 cm from distal end.

TT Silicone s30 Cushion Liner					Durometer : 30 Shore 00 Activity level : (K1-K2)
Liner Size	Length	Circumference 4 cm Above Distal End	Circumference at Proximal End	SKU	Liner Profile
16	35 ± 2 cm	16	20	LNR-SIL-C-[X]-16-30	
18	35 ± 2 cm	18	21.5	LNR-SIL-C-[X]-18-30	
20	35 ± 2 cm	20	23.5	LNR-SIL-C-[X]-20-30	
21	35 ± 2 cm	21	24.5	LNR-SIL-C-[X]-21-30	
22	35 ± 2 cm	22	25.5	LNR-SIL-C-[X]-22-30	
23.5	35 ± 2 cm	23.5	27	LNR-SIL-C-[X]-23.5-30	
25	35 ± 2 cm	25	28.5	LNR-SIL-C-[X]-25-30	
26.5	35 ± 2 cm	26.5	30	LNR-SIL-C-[X]-26.5-30	
28	35 ± 2 cm	28	31.5	LNR-SIL-C-[X]-28-30	
30	35 ± 2 cm	30	33.5	LNR-SIL-C-[X]-30-30	
32	35 ± 2 cm	32	35.5	LNR-SIL-C-[X]-32-30	
34	35 ± 2 cm	34	37.5	LNR-SIL-C-[X]-34-30	
36	35 ± 2 cm	36	39.5	LNR-SIL-C-[X]-36-30	
38	35 ± 2 cm	38	41.5	LNR-SIL-C-[X]-38-30	
40	35 ± 2 cm	40	43.5	LNR-SIL-C-[X]-40-30	
42	35 ± 2 cm	42	45.5	LNR-SIL-C-[X]-42-30	
45	35 ± 2 cm	45	48.5	LNR-SIL-C-[X]-45-30	
TT Silicone s30 Locking Liner					Durometer : 30 Shore 00 Activity level : (K1-K2)
Liner Size	Length	Circumference 4 cm Above Distal End	Circumference at Proximal End	SKU	Liner Profile
16	35 ± 2 cm	16	20	LNR-SIL-L-[X]-16-30	
18	35 ± 2 cm	18	21.5	LNR-SIL-L-[X]-18-30	
20	35 ± 2 cm	20	23.5	LNR-SIL-L-[X]-20-30	
21	35 ± 2 cm	21	24.5	LNR-SIL-L-[X]-21-30	
22	35 ± 2 cm	22	25.5	LNR-SIL-L-[X]-22-30	
23.5	35 ± 2 cm	23.5	27	LNR-SIL-L-[X]-23.5-30	
25	35 ± 2 cm	25	28.5	LNR-SIL-L-[X]-25-30	
26.5	35 ± 2 cm	26.5	30	LNR-SIL-L-[X]-26.5-30	
28	35 ± 2 cm	28	31.5	LNR-SIL-L-[X]-28-30	
30	35 ± 2 cm	30	33.5	LNR-SIL-L-[X]-30-30	
32	35 ± 2 cm	32	35.5	LNR-SIL-L-[X]-32-30	
34	35 ± 2 cm	34	37.5	LNR-SIL-L-[X]-34-30	
36	35 ± 2 cm	36	39.5	LNR-SIL-L-[X]-36-30	
38	35 ± 2 cm	38	41.5	LNR-SIL-L-[X]-38-30	
40	35 ± 2 cm	40	43.5	LNR-SIL-L-[X]-40-30	
42	35 ± 2 cm	42	45.5	LNR-SIL-L-[X]-42-30	
45	35 ± 2 cm	45	48.5	LNR-SIL-L-[X]-45-30	

[X] = Profil thickness in mm (3 or 6)

Example : LNR-SIL-L-**3**-32-30



For a s30 silicone locking liner with 3 mm wall in size 32.

Please see next page for

TT Silicone s40 Liners →



Follow the table to determine the correct size.
Measurements are made 4 cm from distal end.

TT Silicone s40 Cushion Liner					Durometer : 40 Shore 00 Activity level : (K1-K4)
Liner Size	Length	Circumference 4 cm Above Distal End	Circumference at Proximal End	SKU	Liner Profil
16	35 ± 2 cm	16	20	LNR-SIL-C-3-16-40	
18	35 ± 2 cm	18	21.5	LNR-SIL-C-3-18-40	
20	35 ± 2 cm	20	23.5	LNR-SIL-C-3-20-40	
21	35 ± 2 cm	21	24.5	LNR-SIL-C-3-21-40	
22	35 ± 2 cm	22	25.5	LNR-SIL-C-3-22-40	
23.5	35 ± 2 cm	23.5	27	LNR-SIL-C-3-23.5-40	
25	35 ± 2 cm	25	28.5	LNR-SIL-C-3-25-40	
26.5	35 ± 2 cm	26.5	30	LNR-SIL-C-3-26.5-40	
28	35 ± 2 cm	28	31.5	LNR-SIL-C-3-28-40	
30	35 ± 2 cm	30	33.5	LNR-SIL-C-3-30-40	
32	35 ± 2 cm	32	35.5	LNR-SIL-C-3-32-40	
34	35 ± 2 cm	34	37.5	LNR-SIL-C-3-34-40	
36	35 ± 2 cm	36	39.5	LNR-SIL-C-3-36-40	
38	35 ± 2 cm	38	41.5	LNR-SIL-C-3-38-40	
40	35 ± 2 cm	40	43.5	LNR-SIL-C-3-40-40	
42	35 ± 2 cm	42	45.5	LNR-SIL-C-3-42-40	
45	35 ± 2 cm	45	48.5	LNR-SIL-C-3-45-40	
TT Silicone s40 Locking Liner					Durometer : 40 Shore 00 Activity level : (K1-K4)
16	35 ± 2 cm	16	20	LNR-SIL-L-3-16-40	
18	35 ± 2 cm	18	21.5	LNR-SIL-L-3-18-40	
20	35 ± 2 cm	20	23.5	LNR-SIL-L-3-20-40	
21	35 ± 2 cm	21	24.5	LNR-SIL-L-3-21-40	
22	35 ± 2 cm	22	25.5	LNR-SIL-L-3-22-40	
23.5	35 ± 2 cm	23.5	27	LNR-SIL-L-3-23.5-40	
25	35 ± 2 cm	25	28.5	LNR-SIL-L-3-25-40	
26.5	35 ± 2 cm	26.5	30	LNR-SIL-L-3-26.5-40	
28	35 ± 2 cm	28	31.5	LNR-SIL-L-3-28-40	
30	35 ± 2 cm	30	33.5	LNR-SIL-L-3-30-40	
32	35 ± 2 cm	32	35.5	LNR-SIL-L-3-32-40	
34	35 ± 2 cm	34	37.5	LNR-SIL-L-3-34-40	
36	35 ± 2 cm	36	39.5	LNR-SIL-L-3-36-40	
38	35 ± 2 cm	38	41.5	LNR-SIL-L-3-38-40	
40	35 ± 2 cm	40	43.5	LNR-SIL-L-3-40-40	
42	35 ± 2 cm	42	45.5	LNR-SIL-L-3-42-40	
45	35 ± 2 cm	45	48.5	LNR-SIL-L-3-45-40	

Example : LNR-SIL-L-3-22-40
For a s40 silicone locking liner in size 22.

Please see next page for
TT Silicone s30 Liners →



Follow the table to determine the correct size.
Measurements are made 4 cm from distal end.

TF Silicone s40 Locking Liner					Durometer : 40 Shore 00 Activity level : (K1-K3)
Liner Size	Length	Circumference 4 cm Above Distal End	Circumference at Proximal End	SKU	Liner Profil
25	35 ± 2 cm	25	28.5	LNR-TF-SIL-L-2-25-40	
26.5	35 ± 2 cm	26.5	30	LNR-TF-SIL-L-2-26.5-40	
28	35 ± 2 cm	28	31.5	LNR-TF-SIL-L-2-28-40	
30	35 ± 2 cm	30	33.5	LNR-TF-SIL-L-2-30-40	
32	35 ± 2 cm	32	35.5	LNR-TF-SIL-L-2-32-40	
34	35 ± 2 cm	34	37.5	LNR-TF-SIL-L-2-34-40	
36	35 ± 2 cm	36	39.5	LNR-TF-SIL-L-2-36-40	
38	35 ± 2 cm	38	41.5	LNR-TF-SIL-L-2-38-40	
40	35 ± 2 cm	40	43.5	LNR-TF-SIL-L-2-40-40	
42	35 ± 2 cm	42	45.5	LNR-TF-SIL-L-2-42-40	
45	35 ± 2 cm	45	48.5	LNR-TF-SIL-L-2-45-40	
50	35 ± 2 cm	50	58.5	LNR-TF-SIL-L-2-50-40	
55	35 ± 2 cm	55	63	LNR-TF-SIL-L-2-55-40	

Example : LNR-TF-SIL-L-2-38-40
For a TF s40 silicone locking liner in size 38.

Please see next page for
TT Gel liners →



Follow the table to determine the correct size.
Measurements are made 4 cm from distal end.

TT and TF Gel Cushion Liner					Durometer : 30 Shore 00 Activity level : (K1-K3)
Liner Size	Length	Circumference 4 cm Above Distal End	Circumference at Proximal End	SKU	Liner Profil
S	40 ± 2 cm	15-22	24-35	LNR-GEL-C-[X]-S	
M	40 ± 2 cm	19-29	26-40	LNR-GEL-C-[X]-M	
L	40 ± 2 cm	23-37	37-61	LNR-GEL-C-[X]-L	
XL	40 ± 2 cm	27-40	40-64	LNR-GEL-C-[X]-XL	
XXL	40 ± 2 cm	31-53	42-70	LNR-GEL-C-[X]-XXL	
TT and TF Gel Locking Liner					Durometer : 30 Shore 00 Activity level : (K1-K3)
S	40 ± 2 cm	15-22	24-35	LNR-GEL-L-[X]-S	
M	40 ± 2 cm	19-29	26-40	LNR-GEL-L-[X]-M	
L	40 ± 2 cm	23-37	37-61	LNR-GEL-L-[X]-L	
XL	40 ± 2 cm	27-40	40-64	LNR-GEL-L-[X]-XL	
XXL	40 ± 2 cm	31-53	42-70	LNR-GEL-L-[X]-XXL	

[X] = Profil thickness in mm (3, 6 or 9)

Example : LNR-GEL-C-**9**-XXL

For a gel cushion liner with 9 mm wall anterior and 3 mm posterior in size XXL.

Please see next page for

TF Silicone s40 Liners →