

Compatible with pin, lanyard and cushion liners.


To determine the correct size, measure the limb's circumference and length at a specific height.



[Click here](#)

STEP 1


Choose liner type.



PIN

STEP 2


Measure circumference.



Measurement	32-36 cm	36-40 cm	40-44.5 cm	44.5-48 cm	48-52 cm	52-62 cm
SKU	32	38	40	44	48	52

STEP 3

Measure length.



For check sockets you can also measure to the medial trim line.

Measurement	15-20 cm	20-25 cm	25-30 cm	30+ cm
Bladder rows	1	2	3	4
SKU	XSH	SH	LG	XL

Suspension	Circumference	Length	SKU
PIN	32-36 cm	15-20 cm	OVTF-32-XSH-PIN
	36-40 cm		OVTF-38-XSH-PIN
	40-44.5 cm		OVTF-40-XSH-PIN
	44.5-48 cm		OVTF-44-XSH-PIN
	48-52 cm		OVTF-48-XSH-PIN
	52-62 cm		OVTF-52-XSH-PIN
	32-36 cm	20-25 cm	OVTF-32-SH-PIN
	36-40 cm		OVTF-38-SH-PIN
	40-44.5 cm		OVTF-40-SH-PIN
	44.5-48 cm		OVTF-44-SH-PIN
	48-52 cm		OVTF-48-SH-PIN
	52-62 cm		OVTF-52-SH-PIN

Suspension	Circumference	Length	SKU
PIN	32-36 cm	25-30 cm	OVTF-32-LG-PIN
	36-40 cm		OVTF-38-LG-PIN
	40-44.5 cm		OVTF-40-LG-PIN
	44.5-48 cm		OVTF-44-LG-PIN
	48-52 cm		OVTF-48-LG-PIN
	52-62 cm		OVTF-52-LG-PIN
	32-36 cm	30+ cm	OVTF-32-XL-PIN
	36-40 cm		OVTF-38-XL-PIN
	40-44.5 cm		OVTF-40-XL-PIN
	44.5-48 cm		OVTF-44-XL-PIN
	48-52 cm		OVTF-48-XL-PIN
	52-62 cm		OVTF-52-XL-PIN



STEP 1

Choose liner type.

Single seal

1 seal, sealing liners

Overlay will sit under the sealing ring

VAC

Multi seal

4 or 5 seal, sealing liners

Follow instructions for single seal but measure **circumference (C)** just above top seal and **length (L)** from that location to perineum

Overlay will sit above the top sealing ring

STEP 2

Measure circumference.

Circumference

Measure under the reflected sealing ring

Measurement	32-38 cm	38-42 cm	42-47 cm	47-50 cm	50-54 cm	54-62 cm
SKU	32	38	40	44	48	52

STEP 3

Measure length.

Length

Measure under the reflected sealing ring to the perineum

Measurement	20-25 cm	25+ cm
Bladder rows	2	3
SKU	SH	LG

Suspension	Circumference	Length	SKU
VAC	32-38 cm	20-25 cm	OVTF-32-SH-VAC
	38-42 cm		OVTF-38-SH-VAC
	42-47 cm		OVTF-40-SH-VAC
	47-50 cm		OVTF-44-SH-VAC
	50-54 cm		OVTF-48-SH-VAC
	54-62 cm		OVTF-52-SH-VAC

Suspension	Circumference	Length	SKU
VAC	32-38 cm	25+ cm	OVTF-32-LG-VAC
	38-42 cm		OVTF-38-LG-VAC
	42-47 cm		OVTF-40-LG-VAC
	47-50 cm		OVTF-44-LG-VAC
	50-54 cm		OVTF-48-LG-VAC
	54-62 cm		OVTF-52-LG-VAC