

Compatible with distal locking systems, suspension sleeves, cushion liners and lanyards


To determine the most suitable size for the patient, the measurement protocol involves identifying the suspension type, measuring the circumference and length at a specific height, and selecting the appropriate product orientation.



[Click here](#)


STEP 1

Which circumference ?



STEP 2

Which length ?




If patient do not have a socket yet

Length (L)


L = From end of liner to popliteal fossa

(Make sure to measure linear dimensions)




STEP 3

Which side ?




Orientation is based on user preference

Most users choose the same orientation as the side where the Overlay will be worn



Suspension	Orientation	Length (L)	Circumference (C)	Code
Distal locking Suspension sleeve Cushion Lanyards	Right	13.5-20 cm	24-29.5 cm	OV-23-SH-R
			29.5-36 cm	OV-28-SH-R
			36-42 cm	OV-35-SH-R
		+20 cm	24-29.5 cm	OV-23-LG-R
			29.5-36 cm	OV-28-LG-R
			36-42 cm	OV-35-LG-R
	Left	13.5-20 cm	24-29.5 cm	OV-23-SH-L
			29.5-36 cm	OV-28-SH-L
			36-42 cm	OV-35-SH-L
		+20 cm	24-29.5 cm	OV-23-LG-L
			29.5-36 cm	OV-28-LG-L
			36-42 cm	OV-35-LG-L

Compatibility :



Please see next page for
Distal seal →

Compatible with distal seal vacuum systems with or without suspension sleeves


To determine the most suitable size for the patient, the measurement protocol involves identifying the suspension type, measuring the circumference and length at a specific height, and selecting the appropriate product orientation.



[Click here](#)


STEP 1

Which circumference ?




STEP 2

Which length ?




STEP 3

Which side ?



Suspension	Orientation	Length (L)	Circumference (C)	Code
Distal vacuum	Right	13.5-20 cm	24-29.5 cm	OV-23-SH-R
			29.5-36 cm	OV-28-SH-R
			36-42 cm	OV-35-SH-R
		+20 cm	24-29.5 cm	OV-23-LG-R
			29.5-36 cm	OV-28-LG-R
			36-42 cm	OV-35-LG-R
	Left	13.5-20 cm	24-29.5 cm	OV-23-SH-L
			29.5-36 cm	OV-28-SH-L
			36-42 cm	OV-35-SH-L
		+20 cm	24-29.5 cm	OV-23-LG-L
			29.5-36 cm	OV-28-LG-L
			36-42 cm	OV-35-LG-L

Compatibility : Distally positioned seal with or without suspension sleeves



Compatible with :

Össur, Ottobock, Alps, WillowWood liners

And other similar liners

Please see next page for
Other systems →